



DELTA PF CUP: Acetabular Cup
Coppa DELTA PF: Coppa Acetabolare


∅44
LINER SMALL




∅46
LINER SMALL




∅48
LINER SMALL





∅50
LINER MEDIUM



∅52
LINER MEDIUM




∅54
LINER LARGE




DELTA PF CUP: Acetabular Cup
Coppa DELTA PF: Coppa Acetabolare


∅56
LINER LARGE




∅58
LINER LARGE




∅60
LINER LARGE



∅62
LINER LARGE



∅64
LINER LARGE



∅66
LINER LARGE

